

Cell viability, cytotoxicity, proliferation, cell cycle, and apoptosis assays guide

Understand your options when you need
a cell health assay

The analysis of cell viability, cytotoxicity, cell cycle state, cell proliferation, and cell death are critical to most cell-based studies. Use this guide to understand your options when you need to assay cell health.

Choose the right assay to fit your experimental needs:

- Cell viability and cytotoxicity assays based on cell metabolism/enzyme activity

Measure the rate of continuing cellular activities, such as metabolism.

- Cell viability and cytotoxicity assays based on cytolysis/membrane leakage

Test for cell membrane damage, either by measuring the leakage of cellular enzymes or staining with membrane-impermeable dyes.

- Cell proliferation and cell cycle assays

Monitor the growth of a cell population, detect generations of daughter cells, or analyze the cell cycle state of a cell population. You can use these tests to assay cell viability.

- Cell death analysis/apoptosis assays

Measure the markers present in different types of cell death.

Remember: no one method gives a perfect view of cell viability, proliferation, or cell death. It is almost always most effective to combine several different methods.

Cell viability and cytotoxicity assays: cell metabolism and enzyme activity

You can use several methods to assay ongoing cellular metabolism and enzyme activity to quantify cell viability and cytotoxicity.

Methods include

- Dyes reduced by cellular enzymes
- Mitochondrial membrane potential dependent fluorescent dyes
- Cellular esterase cleaved fluorescent dyes
- ATP and ADP assays
- Assays to measure glycolytic flux and oxygen consumption

See below to learn more, or for the deeper analysis of mitochondria, oxidative stress, metabolites, and metabolic activity, see www.abcam.com/metabolismassays.

Dye reduction assays

Tetrazolium cell viability assays rely on cellular dehydrogenases to form a colored formazan product, which is measured by absorbance. Other assays use the reduction of resazurin, by electron acceptance from the mitochondrial respiratory chain, to form the fluorescent resorufin.

Tetrazolium family

Assay	Instrument	Notes	Assay kits
MTT	Plate reader	Original tetrazolium assay; still very popular. Only tetrazolium assay that needs a wash / solubilization step.	ab211091
MTS		Most popular assay. More heavily used than WST-1.	ab197010
WST-1		More sensitive than MTT, XTT or MTS.	ab155902 ab65473 ab65475
Cell Counting Kit-8 / CCK-8 / WST-8			
XTT assay			

Resazurin family

Resazurin is equivalent to the active ingredient of ThermoFisher's alamarBlue®.

Assay	Instrument	Notes	Assay kits
Resazurin	plate reader, microscope, flow cytometer	Fluorometric (Ex/Em 535–560/560–615) or colorimetric. No-wash assay. Fluorescent readout enables multiplexing with other assays.	ab129732




Figure 2. Jurkat cells treated with idarubicin (A) or staurosporin (B) were analyzed with Resazurin assay ab129732.

Mitochondrial membrane potential-dependent dyes

There are several dyes available that accumulate in mitochondria due to the mitochondrial membrane potential and you can use these to identify viable cells. A loss of membrane potential and loss of staining is used to assay for apoptosis.

Assay	Instrument	Notes	Assay kits
TMRE/TMRM	Plate reader, microscope, flow cytometer	Most popular Abcam mitochondrial membrane dye assay. Ex/Em 549/575 nm. Washed out of mitochondria after fixation.	ab113852
JC-1/JC-10		JC-1 (Ex/Em 530/530–570) and JC-10 (Ex/Em 590/520–570) form red aggregates at high concentrations (unaggregated dye is green). Loss of membrane potential causes loss of dye and increased green fluorescence. Washed out after fixation.	JC-1: ab113850 JC-10 (more soluble than JC-1): ab112134 ab112133
Mitotracker Red		Ex/Em 579 /599. Not washed out after fixation.	
Rhodamine 123		Ex/Em 507/529. Washed out after fixation.	
MitoNIR		Plate reader, flow cytometer	Ex/Em 635/660.
MitoOrange	Plate reader, flow cytometer	Ex/Em 540/590.	ab138898 ab138899




Figure 3. Cell staining with TMRE kit ab113852. A: Healthy HeLa cells. B: Healthy Jurkat cells.

Esterase cleavage

Calcein and similar hydrophobic dyes diffuse into cells and are cleaved by intracellular esterases in live cells. The hydrophilic fluorescent product is retained within the cell.

Assay	Instrument	Notes	Assay kits
Calcein AM	Plate reader, microscope, flow cytometer	Ex/Em 495 / 515 nm	
Calcein violet AM	Plate reader, microscope, flow cytometer	Ex/Em 405/460 nm	ab176748
Esterase-cleaved blue	Plate reader	Ex/Em 360/450 nm	ab112120
Esterase-cleaved green	Plate reader, microscope	Ex/Em 490/520 nm	ab112122
Esterase-cleaved near IR	Plate reader	Ex/Em 633/660 nm	ab112123

ATP assays

Most assays use a cell membrane permeabilization agent to release ATP; light is produced using ATP-dependent luciferase. Other ATP assays use the ATP-dependent phosphorylation of glycerol (or other substrates).

Assay	Instrument	Notes	Assay kits
Luminescence ATP assay	Luminometric plate reader	No-wash assay.	ab113849
Luminescence ADP/ATP assay		No-wash assay. After ATP analysis, ADP is converted to ATP for detection.	ab65313
ATP phosphorylation assay	Plate reader	No-wash assay used with cell lysates. Not as sensitive as luminescence assays. Fluorometric (Ex/Em 535/587 nm) is more sensitive than colorimetric.	ab83355

Oxygen consumption and glycolysis assays

The rate of oxygen consumption indicates the level of cellular metabolic activity. Analysis of intracellular oxygen levels and glycolysis activity allow deeper investigation.

Assay	Instrument	Notes	Assay kits
Extracellular oxygen consumption	Plate reader	No-wash assay. Dye signal (Ex/Em 380/650 nm) increases as respiration lowers O ₂ levels. No need for specialized instrument.	ab197243
Intracellular oxygen levels		Dye fluorescence (Ex/Em 340/642) is quenched by intracellular oxygen. No-wash assay.	ab197245
Glycolysis activity		No-wash assay. Lactate production causes extracellular acidification and increased dye fluorescence (Ex/Em 340-380/615 nm).	ab197244




Figure 4. HepG2 cells treated with antimycin A and FCCP and tested with assays for ATP (gray, ab113849), oxygen consumption (white, ab197243) and glycolysis (black, ab197244).

Cell viability and cytotoxicity assays: cytolysis and membrane leakage

You can look at damage to the cell membrane to assay for severe cellular damage or cell death.

Methods include

- Measurement of the activity of enzymes that leak into the extracellular medium
- Membrane impermeable dyes that enter and stain cells upon membrane damage
- Amine-reactive dyes bind weakly to the surface of live cells and create a brighter staining in dead cells

Membrane impermeable dyes are often used, with dyes that stain live cells, in combined live:dead cell assays.

Enzyme leakage

The most popular enzyme leakage cell viability assay is for lactate dehydrogenase.

Assay	Instrument	Notes	Assay kits
LDH/Lactate dehydrogenase	Plate reader	LDH oxidizes lactate and a colored, or fluorescent (Ex/Em 535/587 nm), product is formed.	ab65393 ab197004
AK/Adenylate kinase		AK converts ADP to ATP with detection via luciferase light-generation. AK activity is not as enduring as LDH.	




Figure 5. LDH assay used with staurosporine or cycloheximide HeLa cells, untreated cells, LDH positive control, and lysed cells.

Membrane impermeable dyes

Cell viability assays often use membrane-impermeable fluorescent dyes (mostly DNA stains) that stain cells with damaged cell membranes. Propidium iodide has largely been replaced by DRAQ7™ and 7-AAD for cell viability assays due to its broad emission spectra and tendency to bind to live cells.

Assay	Instrument	Notes	Assay kits
DRAQ7™	Flow cytometer, microscope	Ex/Em 633 & 647/665–800 nm. DNA stain.	ab109202
7-AAD		Ex/Em 488/647 nm. DNA stain.	
Propidium iodide		Ex/Em 536/617 nm. DNA stain. Leaches from cells over time.	ab14083
Ethidium homodimer-1		Ex/Em 528/617. DNA stain.	
Trypan blue	Microscope	Non-fluorescent cell stain. Classic cell viability assay that requires cell counting. Tedious and prone to manual error.	




Figure 6. Jurkat cells treated with staurosporine to induce cell death show DRAQ7™ staining (top half of the plot).

Amine-reactive dyes for live:dead cell assays

Amine-reactive dyes weakly stain viable cells by binding to cell surface amines and strongly stain membrane-compromised cells by reacting with intracellular amines. Dead and live cells can be differentiated by fluorescence level.

Assay	Instrument	Notes	Assay kits
Amine-reactive dyes	Flow cytometer	Ex/Em 410/450 nm. Fixation compatible (applies to all dyes in this table).	ab176738
		Ex/Em 408/512 nm	ab176739
		Ex/Em 398/550 nm	ab176740
		Ex/Em 353/442 nm	ab176741
		Ex/Em 498/521 nm	ab176742
		Ex/Em 547/573 nm	ab176743
		Ex/Em 583/603 nm	ab176744
		Ex/Em 649/660 nm	ab176745

Combined dye live:dead cell assays

Multiple dyes can be combined in a single live:dead cell assay. Examples include the popular Live and dead cell assay (ab115347) with ethidium homodimer to label dead cells and an esterase-cleaved dye for live cells. The alternative, Cell viability assay kit (fluorometric – dual green/red) (ab112121), includes a red DNA staining dye for dead cells and a green esterase-cleaved dye for live cells.





Figure 7. Left: Live:dead cell assay (ab115347) used with live (upper left) and dead (lower right) cells. Right: Etoposide treated cells stained with ab115347. Live cells stain green and dead cells are red.

Cell proliferation and cell cycle assays

These assays monitor the growth rate of a cell population, detect daughter cells in a growing population, or analyze the proportions of cells in different stages of the cell cycle.

Methods include

- DNA staining dyes for cell cycle analysis
- Dye dilution assays
- Incorporation of nucleoside analogs during DNA synthesis

DNA-staining dyes

DNA-staining dyes are commonly used in flow cytometry to measure the DNA content in cell populations and assay for cell cycle state. Propidium iodide is the mostly commonly used dye.

Dye	Instrument	Note	Assay kits
Propidium iodide	flow cytometer	Ex/Em 536/617 nm	ab14083, ab139418
Nuclear Green CCS1		Ex/Em 490/525 nm	ab112116
Nuclear Red CCS1		Ex/Em 490/620 nm	ab112117
DRAQ5™		Ex/Em 633&647/665–800 nm	ab108410
DAPI		Ex/Em 358/461	
Hoechst 33342		Ex/Em 350/461	
Hoechst 33258		Ex/Em 350/461	
7-AAD		Ex/Em 488/647 nm	




Figure 8. Propidium iodide flow cytometry kit (ab139418) used with thymidine (B) and nocodazole (C) treated HeLa cells. Peaks show 2N and 4N DNA content.

Dye dilution assays

The dyes in dye dilution assays are retained within cells over multiple generations. Daughter cells receive half of the dye of parent cells and assays are analyzed on a flow cytometer. Carboxyfluorescein succinimidyl ester (CFSE) is the longest established dye.

Assay	Instrument	Notes	Assay kits
CFSE	flow cytometer	Ex/Em 492/517 nm. Cytotoxic at higher concentrations.	ab113853
CytoLabel Blue	flow cytometer, microscope	Ex/Em 403/454 nm	ab176726
CytoLabel Green		Ex/Em 511/525 nm	ab176735
CytoLabel Red		Ex/Em 628/643 nm	ab176736
CytoLabel Orange		Ex/Em 542/556 nm	ab176737




Figure 9. Flow cytometry analysis of CFSE (ab113853) dilution assay.

Nucleoside analogs

Bromodeoxyuridine (BrdU) and ethynyldeoxyuridine (EdU) assays measure the incorporation of BrdU or EdU into newly synthesized DNA during DNA replication. Unlike BrdU, which is detected using antibodies, EdU can be easily directly labeled, either with a fluorescent dye or biotin for colorimetric or fluorometric detection via streptavidin-HRP. EdU staining is consistent with further antibody staining, unlike the harsher BrdU protocol.

Assay	Instrument	Assay kits
EdU	Microscope, flow cytometry, plate reader	ab219801 ab222421
BrdU	Plate reader, microscope	ELISA: ab126556/ab126572 IHC: ab125306




Figure 10. EdU staining of proliferating HeLa cells. DNA (blue) was stained with Hoechst 33342 (ab145597). Green cells are EdU + Hoechst-positive.

Other methods of assaying proliferation

For the analysis of cell proliferation within tissue samples, or sometimes within cell culture, it is common to use antibodies to stain for the presence of Ki67 or PCNA.

Although little used at high throughput, the classical method of assaying cell proliferation is to use a clonogenic/clonogenicity assay. In this assay, cells are plated out at a low density and then the number of colonies formed is counted.

Senescence assays

The most common marker of senescent cells is the overexpression and accumulation of the endogenous lysosomal beta-galactosidase (SA-beta-gal). Beta-gal activity is detected using a colorimetric or fluorometric substrate.

Assay	Instrument	Assay kits
Beta-gal	microscope, plate reader	ab65351

Cell death/apoptosis assays

You can determine how cells are dying by measuring markers that are activated in different types, and at different stages, of cell death.

Methods include

- Annexin V binding of cell surface phosphatidylserine
- DNA condensation and fragmentation (TUNEL) assays
- Caspase activation and detection assays
- Mitochondrial membrane potential dependent dyes
- Cytochrome C release assays
- Glutathione assays
- Assays for necrosis, anoikis and autophagy

For a more detailed examination of cell death, and assays to analyze it, see our three comprehensive guides to apoptosis, necroptosis, and autophagy [hyperlink to guides].

Annexin V assay

Annexin V binds to phosphatidylserine, which migrates to the outer plasma membrane in apoptosis. Analysis is typically by flow cytometry. Pair Annexin V with a membrane impermeable dye like 7-AAD to distinguish between intact, apoptotic, and necrotic cells (eg see ab214663, ab214484, or ab214485).

Annexin V conjugate	Ex/Em	Assay kits
FITC	495/519	ab14085, ab14082
Cy3	548/561	ab14142, ab14143
Cy5	647/665	ab14150, ab14147
PE	496/576	ab14155, ab14154
PE-Cy5	565/693	ab14159
EGFP	488/530	ab14153, ab14152
Biotin		ab14190, ab14165

For a full list of Annexin V dye conjugates, see www.abcam.com/AnnexinV




Figure 11. Murine platelets were stimulated with CRP in the presence of annexin V-FITC. The annexin V-positive population is indicated (AnV +ve).

DNA condensation/fragmentation

DNA condensation in apoptosis can be measured using DNA stains to visualize condensed nuclei.

DNA fragmentation can be measured using agarose gels. In the TUNEL assay, the 3' ends of DNA fragments are labeled with deoxyuridine either conjugated to a fluorescent dye or biotin.

Assay	Instrument	Assay kits
DNA fragmentation	Gel electrophoresis	ab66090, ab65627, ab66093
TUNEL	Flow cytometry, fluorescence microscope, microscope	ab66108, ab66110, ab206386




Figure 12. TUNEL staining in whole-mount *Hydractinia echinata* using red TUNEL assay kit (ab66110) with DAPI counterstain (blue).

Active caspase detection

Activated caspases can be detected using antibodies with IHC, western blotting, or flow cytometry.

Caspase activity assays either use peptide substrates, which are cleaved by caspases in cell extracts, or similar substrates that bind to activated caspases in live cells. Caspase specificity varies by substrate.

For information on our assays for caspases 1 through 12, formulated either for cell lysates with analysis by plate reader, or for live cells with analysis by flow cytometer, microscope or plate reader, see www.abcam.com/caspaseassays. For more information about caspase substrates, see www.abcam.com/caspsubguide

We also offer assays for cathepsin and calpain activity analysis: cathepsin D (ab65302), cathepsin B (ab65303), cathepsin L (ab65306), and calpains (ab65308).

Mitochondrial membrane potential-dependent dyes

Dyes that accumulate in mitochondria due to the mitochondrial membrane potential are also used in the analysis of apoptosis. For more information, see the earlier section on metabolism-based assays. Apoptotic cells stain more weakly with these dyes due to the loss of membrane potential.

Cytochrome C release

Cytochrome C is released into the cytoplasm following total loss of mitochondrial membrane potential.

Assay	Instrument	Assay kits
Cytochrome C	Western blot, fluorescence microscope	ab110415, ab110417, ab65311

Glutathione assay

Glutathione assays are also used for the analysis of apoptosis.

Assay	Instrument	Assay kits
GSH/GSSG assay	Fluorometric plate reader	ab138881

Necrosis, anoikis, and autophagy

We offer several kits for studying other forms of cell death: necrosis and apoptosis (ab176749 and ab176750), anoikis (ab211153), and autophagy (ab133075 and ab139484).

