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
Sample preparation methods

Cell culture supernatant 細胞培養上清液

- Pipette cell culture media into a centrifuge tube and centrifuge at 1,500 rpm for 10 min at 4°C.
- Immediately aliquot supernatant and store samples at -80°C. Minimize freeze/thaw cycles.

Cell extract 細胞萃取物

- Place tissue culture plates on ice.
- Aspirate medium and gently wash cells once with ice-cold PBS.
- Aspirate PBS and add 0.5 mL complete extraction buffer per 100 mm plate.
- Scrape cells to collect in tilted plate and remove to pre-chilled tube.
- Vortex briefly and incubate on ice for 15–30 min.
- Centrifuge at 13,000 rpm for 10 min at 4°C to pellet insoluble contents.
- Aliquot supernatant (this is the soluble cell extract) to clean, chilled tubes on ice and store samples at -80°C. Minimize freeze/thaw cycles.

Conditioned medium 條件培養基

- Place cells in complete (serum-containing) growth medium and allow cells to proliferate to desired level of confluence.
- Remove growth medium and wash very gently with a few mL of warm PBS. Repeat wash step.
- Remove last PBS wash and gently add serum free growth medium.
- Incubate 1–2 days.
- Pipette medium into a centrifuge tube and centrifuge at 1,500 rpm for 10 min at 4°C.
- Immediately aliquot supernatant and store samples at -80°C. Minimize freeze/thaw cycles.

Milk 牛奶

- Collect samples and centrifuge at 10,000 x g for 2 min at 4°C.
- Aliquot supernatant and store samples at -80°C. Minimize freeze/thaw cycles.

Plasma 血漿

- Collect whole blood into anti-coagulant containing tube or add 0.1 M sodium citrate to 1/10 final volume.
- Centrifuge at 3,000 rpm for 10 min at 4°C.
- Immediately aliquot supernatant (plasma) and store samples at -80°C. Minimize freeze/thaw cycles.

Urine 尿液

- Collect samples and centrifuge at 10,000 x g for 2 min at 4°C.
- Aliquot supernatant and store samples at -80°C. Minimize freeze/thaw cycles.

Saliva 唾液

- Collect samples and centrifuge at 10,000 x g for 2 min at 4°C.
- Aliquot supernatant and store samples at -80°C. Minimize freeze/thaw cycles.

Serum 血清

- Collect whole blood in untreated test tube or, for example, an anti-coagulant free tube.
- Incubate undisturbed at room temperature for 20 min.
- Centrifuge at 3,000 rpm for 10 min at 4°C.
- Immediately aliquot supernatant (serum) and store samples at -80°C. Minimize freeze/thaw cycles.

Tissue extract 組織萃取物

- Dissect the tissue of interest with clean tools, preferably on ice and as quickly as possible to prevent degradation by proteases.
- Place the tissue in round bottom microfuge tubes and immerse in liquid nitrogen to snap freeze. Store samples at -80°C for later use or keep on ice for immediate homogenization.
- For a ~5 mg piece of tissue, add ~300 µL complete extraction buffer (see cell/tissue extraction buffer recipe) to the tube and homogenize with an electric homogenizer.
- Rinse the blade twice using 300 µL complete extraction buffer for each rinse, then maintain constant agitation for 2 h at 4°C (eg place on an orbital shaker in the cold room).
- Centrifuge for 20 min at 13,000 rpm at 4°C. Place on ice, aliquot supernatant (this is the soluble protein extract) to a fresh, chilled tube and store samples at -80°C. Minimize freeze/thaw cycles

Control samples required for ELISA

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We recommend checking the antibody datasheet, which will often provide a suggested positive control.

If no control is suggested, we recommend the following:

- Check to see if there are any Abreviews® for the antibody. Any tissues, cells, or lysates that have been used successfully can be considered a suitable positive control.
- Try looking at the Swiss-Prot or Omnigene database links on the datasheet. These databases will often have a list of tissues that the protein is expressed in. These can also be considered suitable positive controls.
- Check the GeneCards entry for the protein. This will usually provide you with relative levels of expression in various tissues.
- If you still have difficulty finding a suitable control, we recommend doing a quick literature search on PubMed to see which tissues and cells express the protein of interest.

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
Standard in sample matrix (spike control) control

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
Sandwich ELISA protocol

A sandwich ELISA measures antigen between two layers of antibodies (capture and detection antibody). The target antigen must contain at least two antigenic sites capable of binding to antibodies

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- **4- or 5-parameter logistic (4PL or 5PL) curves**

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Scan here to download the full guide!

ELISA troubleshooting tips

Poor standard curve

Cause	Solution
Improper standard solution	Confirm dilutions are made correctly.
Standard improperly reconstituted	Briefly spin vial before opening; inspect for undissolved material after reconstituting.
Standard degraded	Store and handle standard as recommended.
Curve doesn't fit scale	Try plotting using different scales, eg log-log, 5 parameter logistic curve fit.
Pipetting error	Use calibrated pipettes and proper pipetting technique.

No signal

Cause	Solution
Incubation time too short	Incubate samples overnight at 4°C or follow the manufacturer guidelines.
Target present below detection limits of assay	Decrease dilution factor or concentrate samples.
Incompatible sample type	Detection may be reduced or absent in untested sample types. Include a sample that the assay is known to detect a positive control.
Recognition of epitope impeded by adsorption to plate	To enhance detection of a peptide by direct or indirect ELISA, conjugate peptide to a large carrier protein before coating onto the microtiter plate.
Assay buffer compatibility	Ensure assay buffer is compatible with target of interest (eg enzymatic activity retained, protein interactions retained).
Not enough detection reagent	Increase concentration or amount of detection reagent, following manufacturer guidelines.
Sample prepared incorrectly	Ensure proper sample preparation/dilution. Samples may be incompatible with microtiter plate assay format.
Insufficient antibody	Try different concentrations/dilutions of antibody.

Cause	Solution
Incubation temperature too low	Ensure the incubations are carried out at the correct temperature. All reagents including plate should be at room temperature or as recommended by the manufacturer before proceeding.
Incorrect wavelength	Verify the wavelength and read plate again.
Plate washings too vigorous	Check and ensure correct pressure in the automatic wash system. Pipette wash buffer gently if washes are done manually.
Wells dried out	Do not allow wells to become dry once the assay has started. Cover the plate using sealing film or tape for all incubations.
Slow color development of enzymatic reaction	Prepare substrate solution immediately before use. Ensure the stock solution has not expired and is not contaminated. Allow longer incubation.

Large coefficient of variation (CV)

Cause	Solution
Bubbles in wells	Ensure no bubbles are present prior to reading plate.
Wells not washed equally/thoroughly	Check that all ports of the plate washer are unobstructed. Wash wells as recommended.
Incomplete reagent mixing	Ensure all reagents are mixed thoroughly.
Inconsistent pipetting	Use calibrated pipettes and proper technique to ensure accurate pipetting.
Edge effects	Ensure the plate and all reagents are at room temperature.
Inconsistent sample preparation or storage	Ensure consistent sample preparation and optimal sample storage conditions (eg minimize freeze/thaw cycles).

High background

Cause	Solution
Wells are insufficiently washed	Wash wells as per protocol recommendations.
Contaminated wash buffer	Prepare fresh wash buffer.
Too much detection reagent	Ensure the reagent has been diluted properly or decrease the recommended concentration of detection reagent.
Blocking buffer ineffective (eg detection reagent binds blocker; wells not completely blocked)	Try different blocking reagent and/or add blocking reagent to wash buffer.
Salt concentration of incubation/wash buffers	Increasing salt concentrations may reduce non-specific and/or weak off-target interactions.
Waiting too long to read plate after adding stop solution	Read plate immediately after adding stop solution.
Non-specific binding of antibody	Use suitable blocking buffers, eg BSA or 5-10% normal serum - species same as primary antibody if using a directly conjugated detection antibody or same as secondary if using conjugated secondary. Ensure wells are pre-processed to prevent non-specific attachment.
High antibody concentration	Try different dilutions for optimal results.
Substrate incubation carried out in light	Substrate incubations should be carried out in the dark or as recommended by manufacturer.
Precipitate formed in wells upon substrate addition	Increase dilution factor of sample or decrease concentration of substrate.
Dirty plate	Clean the plate bottom.

Low sensitivity

Cause	Solution
Improper storage of ELISA kit	Store all reagents as recommended. Please note that all reagents may not have identical storage requirements.
Not enough target	Concentrate sample or reduce sample dilution.
Inactive detection reagent	Ensure reporter enzyme/fluorophore has the expected activity.
Plate reader settings incorrect	Ensure plate reader is set to read the correct absorbance wavelength or excitation/emission wavelengths for fluorescent detection.
Assay format not sensitive enough	Switch to a more sensitive detection system (eg colorimetric to chemiluminescence/fluorescence). Switch to a more sensitive assay type (eg direct ELISA to sandwich ELISA). Lengthen incubation times or increase temperature.
Target poorly adsorbs to microtiter plate	Covalently link target to microtiter plate.
Not enough substrate	Add more substrate.
Incompatible sample type (eg serum vs cell extract)	Detection may be reduced or absent in untested sample types. Include a sample that the assay is known to detect as a positive control.
Interfering buffers or sample ingredients	Check reagents for any interfering chemicals. For example, sodium azide in antibodies inhibit HRP enzyme and EDTA used as anticoagulant for plasma collection inhibits enzymatic reactions.
Mixing or substituting reagents from different kits	Avoid mixing components from different kits.

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
Everything you need to perform your ELISA experiments

progress happens together
abcam