
# MojoSort™ Dead Cell Removal CITE-Seq Effectively Improved Data Quality

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(A) Sample of fresh isolated PBMCs with 98% live cells, (B) sample of fresh isolated PBMCs spiked to have 40% heat-shocked dead cells, and (C) a sample with spiked-in heat-shocked PBMCs treated with the MojoSort™ Human Dead Cell Removal Kit were stained using TotalSeq™ TBNK lineage antibodies. All cells from samples A, B, and C were collected and pooled for one CITE-Seq experiment. Ridge plots showed a small CD11c peak in all samples with the exception of a middle peak from sample B. This unexpected peak (circled in red) can be eliminated after dead cell removal. (E) Cells with higher percent mitochondrial reads were considered dead and dying cells. These apoptotic and dead cells are abundant in sample B, and can be eliminated by the use of Dead Cell Removal Kit as in sample C. (F) UMAP based on TBNK lineage protein marker expression showed that common PBMC cell types can be detected from all three samples, with the exception of a dead cell cluster in the center that consists of mostly cells from the artial dead cell sample B. (G) Ridge plots showed that sample C preserved all the common lineages in PBMCs after dead cell removal, and the expression levels of the lineage markers were mostly identical when comparing to samples A and C.