

糖尿病與心血管疾病研究試劑盒

Human

Product Name	Assay Range	Sensitivity	Cat No.
Human Proinsulin ELISA	3.3-132 pmol/L	0.5 pmol/L	10-1118-01
Human C-peptide ELISA	100-4000 pmol/L	15 pmol/L	10-1136-01
Human C-peptide ELISA (Ultrasensitive)	5-230 pmol/L	1.5 pmol/L	10-1141-01
Iso-Insulin ELISA	3-100 mU/L	1 mU/L	10-1128-01
 Human Insulin ELISA (Ultrasensitive)	0.15-20 mU/L	0.07 mU/L	10-1132-01
Human/Mouse/Rat Glucagon ELISA	1.5-130 pmol/L	1.5 pmol/L	10-1271-01
Human Glicentin ELISA	3-300 pmol/L	3 pmol/L	10-1273-01
Total GIP NL-ELISA	2.7-1000 pmol/L	1.62 pmol/L	10-1258-01
Mercodia Total GLP-1 NL-ELISA	0.9 - 940 pmol/L	0.65 pmol/L	10-1278-01
Lispro NL-ELISA	1-500 mU/L	0.36 mU/L	10-1291-01

Rat & Mouse

Product Name	Assay Range	Sensitivity	Cat No.
 Mouse Insulin ELISA	0.2-6.5 µg/L	0.2µg/L	10-1247-01
Rat/Mouse Proinsulin ELISA	3.5-200 pmol/L	3 pmol/L	10-1232-01
Rat Insulin ELISA	0.15-5.5 µg/L	0.15µg/L	10-1250-01
Rat Insulin ELISA (High Range)	3-150 µg/L	1.5µg/L	10-1145-01
Rat Insulin ELISA (Ultrasensitive)	0.02-1 µg/L	0.020µg/L	10-1251-01
Rat/Mouse Glucagon ELISA	2-180 pmol/L	1.5 pmol/L	10-1281-01
Rat C-peptide ELISA	100-4000 pmol/L	100 pmol/L	10-1172-01



Larger animals

Product Name	Assay Range	Sensitivity	Cat No.
Porcine C-peptide ELISA	20-1200 pmol/L	10 pmol/L	10-1256-01
Porcine Insulin ELISA	0.02-1.5 µg/L	0.01µg/L	10-1200-01
Bovine Insulin ELISA	0.05-3.0 µg/L	0.025µg/L	10-1201-01
Canine Insulin ELISA	0.02-1.5 µg/L	0.01µg/L	10-1203-01
Equine Insulin ELISA	0.02-1.5 µg/L	0.01µg/L	10-1205-01
Feline Insulin ELISA	10-700 ng/L	9.2 ng/L	10-1233-01
Ovine Insulin ELISA	0.05-3.0 µg/L	0.025µg/L	10-1202-01

Cardiovascular diseases

Product Name	Assay Range	Sensitivity	Cat No.
Human MPO ELISA	3.5-200 µg/L	3 µg/L	10-1176-01
Human Oxidized LDL ELISA	1.4-22.5 mU/L	< 1 mU/L	10-1143-01
Human Lp(a) ELISA	0.3-5.0 U/L	0.05 U/L	10-1106-01

昶安科技 EXCEL BIOMEDICAL

www.exbio.com.tw

0800-0123-10

service@exbio.com.tw

