


Mitochondrion

粒 腺 體



MitoSpy™ 粒線體定位探針是可滲透細胞的螢光化學試劑，用於標記活細胞中的粒線體。
MitoSpy™ 根據其膜電位定位到粒線體，可用於指示細胞健康狀況以及定位。
它還可以被固定和透化以用於進一步的基於抗體的檢測。

品名	反應物種	應用	貨號
JC-10 Mitochondrial Membrane Potential Kit	Human, Mouse, Rat, All Species	ICC (live imaging), FC	421902
MitoSpy™ Green FM	Human, Mouse, Rat, All Species	ICC	424805
MitoSpy™ Orange CMTMRos	Human, Mouse, Rat, All Species	ICC, FC	424803
MitoSpy™ Red CMXRos	Human, Mouse, Rat, All Species	ICC, FC	424801
MitoSpy™ NIR DiIC1 (5)	Human, Mouse, Rat, All Species	ICC (live imaging), FC	424807