

RNA Isolation Kit

BioVendor RNA isolation kit 可用於純化全RNA，
適用於全血、血漿、血清、培養細胞和組織等。



Total RNA purification

small RNA (<200 nt) &
large RNA (>200 nt) fractions



Fast isolation

approximately 45 minutes



Excellent yield

High recovery and purity



User-friendly

No phenol/chloroform extraction



產品規格

- Silica-membrane technology
- Mini spin columns
- Fragment size: ≥ 18 nt
- Binding capacity 200 μ g
- Elution volume 20 - 100 μ l



樣本需求

- 300 μ L plasma or serum
- <107 cultured cells
< 30 mg human/animal tissue
- <50 mg plant tissue



產品表

品名	規格	貨號
Cell/Exosome miRNA isolation kit	50 rxn	RMR001R
Tissue miRNA isolation kit	50 rxn	RMR004R
Serum/plasma miRNA isolation kit	50 rxn	RMR002R
FFPE tissue miRNA isolation kit	50 rxn	RMR003R

