

Steroid Hormone Assay

類固醇激素試劑盒

Cayman 為脂質類化學物質的專家，推出多種類固醇激素試劑盒，可用於免疫、代謝、發炎、再生、壓力性等試驗，多年來深受科學家喜愛。



產品優勢



高靈敏度



適用於各個物種



引用文獻高

Androgens

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
582701	Testosterone ELISA Kit	3.9-500 pg/ml	6 pg/ml
582751	11-keto Testosterone ELISA Kit	0.78-100 pg/ml	1.3 pg/ml
501740	Testosterone Glucuronide ELISA Kit	4.9-3,000 pg/ml	25 pg/ml

Estrogens

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
501890	Estradiol ELISA Kit	0.61-10,000 pg/ml	20 pg/ml
582281	Estriol ELISA Kit	0.9-2,000 pg/ml	4 pg/ml (4°C for 18 hrs), 12 pg/ml (RT for 2 hrs)
501290	Estrone-3-Glucuronide ELISA Kit	2.7-6,000 pg/ml	23.2 pg/ml (overnight), 36.1 pg/ml (rapid)
582261	2-Methoxyestradiol ELISA Kit	13.1-8,000 pg/ml	40 pg/ml

Progestogens

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
501300	Pregnanediol-3-Glucuronide (PDG) ELISA Kit	0.4-50 ng/ml	0.94 ng/ml
582601	Progesterone ELISA Kit	7.8-1,000 pg/ml	10 pg/ml

Glucocorticoids

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
501320	Corticosterone ELISA Kit	8.2-5,000 pg/ml	30 pg/ml
500360	Cortisol ELISA Kit	6.6-4,000 pg/ml	35 pg/ml
500370	Cortisol Express ELISA Kit	39.1-5,000 pg/ml	110 pg/ml
501420	11-Oxoetiocholanolone ELISA Kit	0.061-1,000 ng/ml	0.2 ng/ml

Mineralocorticoids

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
501090	Aldosterone ELISA Kit	15.6-2,000 pg/ml	30 pg/ml

Vitamin D

Item No.	Product Name	Assay Range	Sensitivity (80% B/BQ)
501050	Vitamin D ELISA Kit	0.19-25 ng/ml	0.5 ng/ml



Steroid Hormone Receptor

類固醇激素受體

類固醇激素通過結合併激活特定的核受體發揮其大部分作用，這些核受體是細胞內配體誘導的轉錄因子。

Cayman 提供基於細胞的檢測試劑盒來測量多個物種的類固醇激素受體活性。



Androgen Receptor

Item No.	Product Name
15756	Human Androgen Receptor Reporter Assay System
28338	Rat Androgen Receptor Reporter Assay System
32578	Zebrafish Androgen Receptor Reporter Assay System

Estrogen Receptors (ERs)

ER α

Item No.	Product Name
15737	Human Estrogen Receptor Alpha Reporter Assay System
36709	Rat Estrogen Receptor Alpha Reporter Assay System
30256	Zebrafish Estrogen Receptor Alpha Reporter Assay System

ER β

Item No.	Product Name
15738	Human Estrogen Receptor Beta Reporter Assay System
36710	Rat Estrogen Receptor Beta Reporter Assay System

ER α & ER β

Item No.	Product Name
15739	Human Estrogen Receptors Reporter Assay Panel

Progesterone Receptor

Item No.	Product Name
15759	Human Progesterone Receptor Reporter Assay System
32581	Rat Progesterone Receptor Reporter Assay System

Glucocorticoid Receptor

Item No.	Product Name
15733	Human Glucocorticoid Receptor Reporter Assay System
21400	Mouse Glucocorticoid Receptor Reporter Assay System
33510	Rat Glucocorticoid Receptor Reporter Assay System
32579	Zebrafish Glucocorticoid Receptor Reporter Assay System

Mineralocorticoid Receptor

Item No.	Product Name
15740	Human Mineralocorticoid Receptor Reporter Assay System

Vitamin D Receptor

Item No.	Product Name
15742	Human Vitamin D Receptor Reporter Assay System



類固醇激素
常見問題



Cayman
年度熱銷產品



Cayman
氧化壓力試劑盒

